

Key TLR Ligands

Tomorrow's Reagents Manufactured Today®

International Edition

For TLR2 & TLR6/2

MALP-2

ALX-162-027-C050
ALX-162-034-C500

50 µg
500 µg

BULK

Synthetic. Originally isolated from *Mycoplasma fermentans*. This MALP-2 corresponds to the originally isolated isomer, which expresses potent endotoxin-like activity and approaches in certain experimental systems the toxicity of LPS. MALP-2 signalling, unlike that of LPS, is induced via TLR2 and TLR6 signalling.

LIT: Purification and partial biochemical characterization of a *Mycoplasma fermentans*-derived substance that activates macrophages to release nitric oxide, tumor necrosis factor, and interleukin-6: PF. Muhlradt and M. Frisch; *Infect. Immun.* **62**, 3801 (1994) • For a comprehensive bibliography please visit our website.

For more Details see Page 3.

For TLR3

Polyinosinic-polycytidylic acid . K

[poly(I:C) . K (TLRgrade™) (synthetic)]

ALX-746-021-M005

5 mg

Specific ligand for TLR3 [1, 2] and MDA5/Helicard [3].

LIT: [1] The dsRNA binding site of human Toll-like receptor 3: J.K. Bell, et al; *PNAS* **103**, 8792 (2006) • [2] Subcellular localization of Toll-like receptor 3 in human dendritic cells: M. Matsumoto, et al; *J. Immunol.* **171**, 3154 (2003) • [3] Differential roles of MDA5 and RIG-I helicases in the recognition of RNA viruses: H. Kato, et al; *Nature* **441**, 101 (2006)

For TLR5

Flagellin (high purity)

ALX-522-058-C010

10 µg

Isolated from *Salmonella typhimurium* strain 14028. **SPECIFICITY:** Binds to human and mouse TLR5. **BIOLOGICAL ACTIVITY:** Activation of TLR5 in human epithelial cell assays based on NF-κB luciferase fusions.

LIT: Flagellin stimulation of intestinal epithelial cells triggers CCL20-mediated migration of dendritic cells: F. Sierro, et al; *PNAS* **98**, 13722 (2001) • Pathophysiological role of Toll-like receptor 5 engagement by bacterial flagellin in colonic inflammation: S.H. Rhee, et al; *PNAS* **102**, 13610 (2005)

For TLR7 & 8

R-848

[S 28463; Resiquimod]

ALX-420-038-M005

5 mg

ALX-420-038-M025

25 mg

Selective ligand for TLR7 in mouse and for TLR7 and TLR8 in human. Potent antitumor and antiviral compound.

LIT: The immune response modifier resiquimod mimics CD40-induced B cell activation: G.A. Bishop, et al; *Cell. Immunol.* **208**, 9 (2001) • Small anti-viral compounds activate immune cells via the TLR7 MyD88-dependent signaling pathway: H. Hemmi, et al; *Nat. Immunol.* **3**, 196 (2002) • For a comprehensive bibliography please visit our website.

Also available:

TLR4 Ligands: Ultra-pure TLRgrade™ LPS. For more information see Pages 4-5.

TLR9 Ligands: Ultra-pure TLRgrade™ ODNs. For more information see Pages 6-7.

For a complete Overview on TLR Ligands see Pages 2-3.

highlight

Toll-like Receptors & Their Ligands

	Ligands	Origin
TLR1 / TLR2	TLR1 complexed with TLR2	
ALX-165-066/-069	Lipoproteins/triacylated lipopeptides: Pam₃CSK₄ , JBT3002, OspA	Bacteria, Mycobacteria
	Soluble Lipoproteins	<i>Neisseria meningitidis</i>
TLR2		
	Bacterial Lipoproteins (BLPs)	Bacteria
	Lipoarabinomannan (LAM)	Mycobacteria
ALX-162-027	MALP-2 (Mycoplasmal Macrophage-activating Lipopeptide-2)	Mycoplasma
	Glycosylphosphatidylinositol (GPI)	<i>Trypanosoma cruzi</i>
	Glycolipids	<i>Treponema maltophilum</i>
	Porins	<i>Neisseria sp.</i>
TLR2/Dectin-1	TLR2 interacts with Dectin-1	
	Zymosan	Fungi
TLR3		
	Double-stranded RNA (dsRNA)	Viral
ALX-746-021	Polyinosine-polycytidylic Acid (poly(I:C))	Synthetic
	mRNA	Host
	tRNA	Host/fungi?
TLR4 / CD14	TLR4/MD-2 complexed with CD14 and/or LPS-binding Protein [LBP]	
ALX-581-007 to ALX-581-020 & ALX-581-150	S-Lipopolysaccharides (LPS) (smooth) wild-type (wt) LPS (contains repeated O-polysaccharide units)	Gram-negative bacteria (e.g. <i>E. coli</i> , <i>Salmonella</i>)
	R-Lipopolysaccharides (LPS) (rough) mutant (Ra, Rb: extended core-polysaccharide) LPS (S-LPS-like)	Gram-negative bacteria (e.g. <i>E. coli</i> , <i>Salmonella</i>)
	R-Lipopolysaccharides (LPS) (rough) mutant (Rc, Rd1, Rd2, Re: short core-polysaccharide) LPS (Lipid A-like)	Gram-negative bacteria (e.g. <i>E. coli</i> , <i>Salmonella</i>)
	Flavolipin	<i>Flavobacterium meningosepticum</i>
ALX-351-001	Paclitaxel (Taxol®) (recognizes mouse TLR4)	Plant
	Acyclic Lipid A-like Analog (R-112022)	Synthetic
	Type III Repeat Extra Domain A (EDA)(*)	Host
	LMW Oligosaccharides of Hyaluronic Acid (sHA)	Host
	Polysaccharide Fragments of Heparan Sulfate(*)	Host (only mouse tested)
	Fibrinogen(*)	Host
	Fusion Protein of RSV(*)	Respiratory syncytial virus
	Envelope Proteins of MMTV(*)	Mouse mammary tumor virus
	Glycoinositolphospholipids (GIPLs)	<i>Trypanosoma cruzi</i>
	Heat Shock Proteins (HSPs)(*)	TLR2 vs TLR4
TLR4	TLR4/MD-2 in the absence of CD14 and/or LPS-binding Protein [LBP]	
see above	R-Lipopolysaccharides (LPS) (rough) mutant (Rc, Rd1, Rd2, Re: short core-polysaccharide) LPS (Lipid A-like)	Gram-negative bacteria (e.g. <i>E. coli</i> , <i>Salmonella</i>)
ALX-581-200 to -203	Lipid A (lacks core-polysaccharide: active component of LPS)	Gram-negative bacteria (e.g. <i>E. coli</i> , <i>Salmonella</i>)
TLR5		
ALX-522-058	Flagellin	Gram-negative bacteria
TLR6 / TLR2	TLR6 complexed with TLR2	
ALX-162-027	MALP-2 (Mycoplasmal Macrophage-activating Lipopeptide-2) / Diacylated Macrophage-activating Lipopeptide-2	Mycoplasma
	Diacylated Lipopeptide FSL-1	Part of the lipoprotein LP44 of <i>Mycoplasma salivarium</i>
	Diacylated Lipopeptide Pam ₂ CSK ₄	Synthetic
	Soluble Tuberculosis Factor (STF)	Mycobacteria
TLR7		
ALX-420-039/-040	Imiquimod / Gardiquimod	Synthetic
ALX-420-038	R-848 (Resiquimod)	Synthetic
	S-27610	Synthetic
ALX-480-097	Loxoribine / TOG	Synthetic
	3M-13	Synthetic
	Bropirimine	Synthetic
	Single-stranded RNA (ssRNA) / Polyuridylic acid . K (Poly(U) . K)	Non-viral & viral ssRNA
	U1snRNA	Host
	siRNA	Synthetic
TLR8		
	3M-2	Synthetic
ALX-420-038	R-848 (Resiquimod)	Synthetic
	Single-stranded RNA (ssRNA) / Polyuridylic acid . K (Poly(U) . K)	Non-viral & viral ssRNA
TLR9	TLR9 (intracellular and extracellular)	
	Unmethylated CpG DNA	Bacteria / Protozoa/ viral
ALX-746-001 to -006	CpG ODNs	Synthetic
ALX-746-024 to -026	AT-ODNs	Synthetic
	Hemozoin	Protozoa (<i>Malaria piment</i>)
TLR10		
	Unknown	
TLR11		
ALX-522-093	Profilin	Protozoa (<i>Toxoplasma gondii</i>)
	Uropathogenic <i>E. coli</i> (UPECs)	Bacteria

(*) See Technical Note on Page 4

Key Toll-like Receptor Ligands

TLR1/2

Pam₃Cys-Ser-(Lys)₄ . 3HCl

[Pam₃CSK₄ . 3HCl]


ALX-165-066-M002

2 mg

Selective agonist of TLR1 complexed with TLR2. Cell permeable, water soluble synthetic cationic lipohexapeptide analog of the immunologically active N-terminal portion of bacterial lipoprotein that potently activates monocytes and macrophages.

LIT: Lipopeptide derivatives of bacterial lipoprotein constitute potent immune adjuvants combined with or covalently coupled to antigen or hapten: A. Reitermann, et al.; Biol. Chem. Hoppe Seyler 370, 343 (1989)


■ For a comprehensive bibliography please visit our website.



Pam₃Cys-Ser-(Lys)₄ (Aca-Aca-Biotin) . 2TFA

ALX-165-069-M001

Selective, biotin-labelled ligand of TLR1 complexed with TLR2.



TLR2 & TLR6/2

MALP-2

ALX-162-027-C050

50 µg

ALX-162-034-C500

500 µg

BULK

Synthetic. MALP-2 was originally isolated from *Mycoplasma fermentans*. This MALP-2 corresponds to the originally isolated isomer, which expresses potent endotoxin-like activity and approaches in certain experimental systems the toxicity of LPS. MALP-2 signalling, unlike that of LPS, is not transduced via TLR4, but is induced via TLR2 and TLR6 signalling. For more information about MALP-2 (incl. stereochemistry) see www.malp-research.de. **SPECIFIC ACTIVITY:** ~2x10⁸units/mg. One unit is defined as the dilution giving half maximal release of nitric oxide from C3H/HeJ mouse peritoneal exudate cells in the standard assay.

LIT: Purification and partial biochemical characterization of a *Mycoplasma fermentans*-derived substance that activates macrophages to release nitric oxide, tumor necrosis factor, and interleukin-6: P.F. Muhlradt and M. Frisch; Infect. Immun. 62, 3801 (1994) ■ For a comprehensive bibliography please visit our website.

TLR3 –

Specific Ligand for TLR3

Polyinosinic-polycytidylic acid . K

[poly(I:C) . K (TLRgrade™) (synthetic)]

ALX-746-021-M005

5 mg

Specific ligand for TLR3 [1, 2] and MDA5/Helicard [3].

LIT: [1] The dsRNA binding site of human Toll-like receptor 3: J.K. Bell, et al.; PNAS 103, 8792 (2006) ■ [2] Subcellular localization of Toll-like receptor 3 in human dendritic cells: M. Matsumoto, et al.; J. Immunol. 171, 3154 (2003) ■ [3] Differential roles of MDA5 and RIG-I helicases in the recognition of RNA viruses: H. Kato, et al.; Nature 441, 101 (2006)

TLR4

Paclitaxel

[Taxol®]

ALX-351-001-M001

1 mg


ALX-351-001-M005

5 mg

ALX-351-001-M025

25 mg

Isolated from *Taxus brevifolia*. LPS mimetic in mouse but not human involving the TLR4 signalling pathway.



TLR5

Flagellin (high purity)

ALX-522-058-C010

10 µg

Isolated from *Salmonella typhimurium* strain 14028. **SPECIFICITY:** Binds to human and mouse TLR5. **BIOLOGICAL ACTIVITY:** Activation of TLR5 in human epithelial cell assays based on NF-κB luciferase fusions.

LIT: Flagellin stimulation of intestinal epithelial cells triggers CCL20-mediated migration of dendritic cells: F. Sierro, et al.; PNAS 98, 13722 (2001) ■ Pathophysiological role of Toll-like receptor 5 engagement by bacterial flagellin in colonic inflammation: S.H. Rhee, et al.; PNAS 102, 13610 (2005)

TLR7/8

Gardiquimod

ALX-420-040-M025

25 mg

ALX-420-040-M100

100 mg

Selective ligand for human or mouse TLR7. Induces the activation of NF-κB in HEK 293 cells expressing TLR7. At high concentrations (3µg/ml) slightly activates TLR8. More active than imiquimod (Prod. No. ALX-420-039).

Imiquimod

[R-837]

ALX-420-039-M100


100 mg

ALX-420-039-M250

250 mg

Topical immune response modifier that inhibits angiogenesis. Up-regulates IL-18 and down-regulates MMP-9 through recognition of TLR7 and subsequent activation of MyD88-dependent pathway.

LIT: Imiquimod applied topically: a novel immune response modifier and new class of drug: R.L. Miller, et al.; Int. J. Immunopharmacol. 21, 1 (1999) ■ For a comprehensive bibliography please visit our website.



Loxoribine

[7-Allyl-8-oxoguanosine]

ALX-480-097-M025


25 mg

ALX-480-097-M100

100 mg

Selective ligand for TLR7 via IFN. Activates natural killer cells and primes cytolytic precursor cells for activation by IL-2.

LIT: Loxoribine (7-allyl-8-oxoguanosine) activates natural killer cells and primes cytolytic precursor cells for activation by IL-2: B.L. Pope, et al.; J. Immunol. 151, 3007 (1993) ■ For a comprehensive bibliography please visit our website.



R-848

[S 28463; Resiquimod]

ALX-420-038-M005

5 mg

ALX-420-038-M025


25 mg

Selective ligand for TLR7 in mouse and for TLR7 and TLR8 in human. Potent antitumor and antiviral compound.

LIT: The immune response modifier resiquimod mimics CD40-induced B cell activation: G.A. Bishop, et al.; Cell. Immunol. 208, 9 (2001)

■ Small anti-viral compounds activate immune cells via the TLR7 MyD88-dependent signaling pathway: H. Hemmi, et al.; Nat. Immunol. 3, 196 (2002)

■ For a comprehensive bibliography please visit our website.



TLR11

Profilin (*Toxoplasma gondii*) (rec.)

ALX-522-093-C010

10 µg

Produced in *E. coli*. Full length profilin (aa 1-163) from *Toxoplasma gondii* is fused to a N-terminal tag. **SPECIFICITY:** Binds to mouse TLR11.

TLR Ligand Sets

TLR Ligands Set I

APO-54N-018-KI01

1 Set

KIT/SET CONTAINS: Sterile, ready-to-use solutions of ligands to TLRs 1 to 9: 25µg of Pam₃Cys-Ser-(Lys)₄ . 3HCl (TLR1/2) / 1 mg Polyinosinic-polycytidylic acid . potassium salt (TLR3) / 20µg LPS from *E. coli*, Serotype R515 (Re) (TLR4) / 2µg Flagellin (purified) (TLR5) / 2µg MALP-2 (TLR6/2) / 100µg Polyuridylic acid . potassium salt (TLR7/8) / 100µg CpG ODN 2395 (TLR9). For activation of TLRs 1 to 9.

LIT: A soluble form of lymphocyte activation gene-3 (IMP321) induces activation of a large range of human effector cytotoxic cells: C. Brignone, et al.; J. Immunol. 179, 4202 (2007)

TLR (human) Ligands Set II

APO-54N-030-KI01

1 Set

KIT/SET CONTAINS: Sterile, ready-to-use solutions of ligands to TLRs: 25µg of Pam₃Cys-Ser-(Lys)₄ . 3HCl (TLR1/2) / 1 mg Polyinosinic-polycytidylic acid . potassium salt (TLR3) / 50µg LPS from *E. coli*, Serotype R515 (Re) (TLR4) / 10µg Flagellin (purified) (TLR5) / 2µg MALP-2 (TLR6/2) / 100µg Imiquimod (TLR7) / 100µg R848 (TLR7/8) / 100µg CpG ODN 2216 (Type A) (TLR9). For activation of human TLRs.

TLR (mouse) Ligands Set III

APO-54N-031-KI01

1 Set

KIT/SET CONTAINS: Sterile, ready-to-use solutions of ligands to TLRs: 25µg of Pam₃Cys-Ser-(Lys)₄ . 3HCl (TLR1/2) / 1 mg Polyinosinic-polycytidylic acid . potassium salt (TLR3) / 50µg LPS from *E. coli*, Serotype R515 (Re) (TLR4) / 10µg Flagellin (purified) (TLR5) / 2µg MALP-2 (TLR6/2) / 100µg Imiquimod (TLR7) / 100µg CpG ODN 1585 (Type A) (TLR9) / 10µg Profilin (from *Toxoplasma gondii*) (TLR11). For activation of mouse TLRs 1 to 11.

TLRgrade™ – LPS & Lipid A Reagents

Sterile, ready-to-use liquid formulation – No hazardous handling! No further purification required! Potent and selective activators of TLR4 (10–1000 ng/ml). Ultrapure (≥99.9%) – Concentration: 1 mg/ml

Prod. No.	Product	Size
ALX-581-015-L002	LPS from <i>Salmonella minnesota</i> R345 (Rb) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-017-L002	LPS from <i>Salmonella minnesota</i> R5 (Rc) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-018-L002	LPS from <i>Salmonella minnesota</i> R7 (Rd) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-008-L002	LPS from <i>Salmonella minnesota</i> R595 (Re) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-010-L002	LPS from <i>E. coli</i> , Serotype EH100 (Ra) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-014-L002	LPS from <i>E. coli</i> , Serotype J5 (Rc) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-007-L002	LPS from <i>E. coli</i> , Serotype R515 (Re) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-012-L002	LPS from <i>E. coli</i> , Serotype O111:B4 (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-013-L002	LPS from <i>E. coli</i> , Serotype O55:B5 (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-009-L002	LPS from <i>Salmonella abortus equi</i> S-form (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-150-R500	NEW LPS from <i>Salmonella abortus equi</i> S-form (TLRgrade™) (Ready-to-Use) (Biotin)	500 µl
ALX-581-011-L002	LPS from <i>Salmonella typhimurium</i> S-form (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-019-L002	NEW LPS from <i>Salmonella enteritidis</i> S-form (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-020-L002	NEW LPS from <i>Salmonella minnesota</i> S-form (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-200-L002	Lipid A from <i>E. coli</i> , Serotype R515 (Re) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-201-L002	Lipid A from <i>Salmonella minnesota</i> R595 (Re) (TLRgrade™) (Ready-to-Use)	2 ml
ALX-581-203-L001	NEW Monophosphoryl Lipid A [MPL-A] from <i>E. coli</i> , Serotype R515 (Re) (TLRgrade™) (Ready-to-Use)	1 ml
ALX-581-202-L001	Monophosphoryl Lipid A [MPL-A] from <i>Salmonella minnesota</i> R595 (Re) (TLRgrade™) (Ready-to-Use)	1 ml

All the TLRgrade™ LPS and Lipid A preparations are specific activators of Toll-like receptor (TLR) 4 and **do not activate TLR2 or other TLRs** as determined with splenocytes and macrophages from TLR4 deficient mice (see Figure 1).




FIGURE 1: Activation of macrophages from TLR4 wild type compared to TLR4 deficient mice by TLRgrade™ LPS and TLRgrade™ Lipid A from ALEXIS® Biochemicals. Lipid A or LPS concentrations, which induced maximal activation of TLR4 wild type mouse macrophages, were also applied to TLR4 deficient mouse macrophages (10⁵) as follows: 80ng *Salmonella* LPS (*S. minnesota* R595 Re), 80ng *E. coli* LPS (R515 Re), 400ng *E. coli* Lipid A (R515 Re) and 400ng *Salmonella* LPS (*S. abortus equi* S-form). 10 units of IL-6 correspond to the detection limit of the IL-6 ELISA.

Technical Note

Contaminated TLR Ligands

The delicate and specific recognition of different PAMPs by TLRs revealed that only the purest ligands, free of any other immunostimulatory contamination, allow to successfully elucidate the role of each TLR. **While LPS was thought to not only activate TLR4 but also TLR2, repurification of commercial preparations of both *E. coli* and *Salmonella minnesota* showed that LPS no longer induces cellular activation through TLR2 [1-5].** On the other hand highly purified HSP60 [6-7] and HSP70 [8] do not stimulate TLR4 as previously reported [9]. Furthermore it has been shown that purified peptidoglycans activate Nod1 and does not involve TLR2 or TLR4 [10, 11]. Even synthetic CpG ODNs show different activation of certain immune cell subsets when highly purified (TLRgrade™ CpG ODNs, see next chapter).

LIT: [1] Lipopolysaccharides (LPS) of oral black-pigmented bacteria induce tumor necrosis factor production by LPS-refractory C3H/HeJ macrophages in a way different from that of *Salmonella* LPS: T. Kirikae, et al.; *Infect. Immun.* **67**, 1736 (1999) • [2] Repurification of lipopolysaccharide eliminates signalling through both human and murine toll-like receptor 2: M. Hirschfeld, et al.; *J. Immunol.* **165**, 618 (2000) • [3] Toll-like receptor 4, but not toll-like receptor 2, is a signalling receptor for *Escherichia* and *Salmonella* lipopolysaccharides: R.I. Tapping, et al.; *J. Immunol.* **165**, 5780 (2000) • [4] Two lipoproteins extracted from *Escherichia coli* K-12 LCD25 lipopolysaccharide are the major components responsible for Toll-like receptor 2-mediated signalling: H.K. Lee, et al.; *J. Immunol.* **168**, 4012 (2002) • [5] Murein lipoprotein, peptidoglycan-associated lipoprotein, and outer membrane protein A are present in purified rough and smooth lipopolysaccharides: J. Hellman, et al.; *J. Infect. Dis.* **188**, 286 (2003) • [6] Recombinant human heat shock protein 60 does not induce the release of tumor necrosis factor alpha from murine macrophages: B. Gao & M.F. Tsan; *J. Biol. Chem.* **278**, 22523 (2003) • [7] Endotoxin-free heat-shock protein 70 fails to induce APC activation: H. Bausinger, et al.; *Eur. J. Immunol.* **32**, 3708 (2002) • [8] Endotoxin contamination in recombinant human heat shock protein 70 (Hsp70) preparation is responsible for the induction of tumor necrosis factor alpha release by murine macrophages: B. Gao & M.F. Tsan; *J. Biol. Chem.* **278**, 174 (2003) • [9] Interaction of TLR2 and TLR4 ligands with the N-terminal domain of Gp96 amplifies innate and adaptive immune responses: T. Warger, et al.; *J. Biol. Chem.* **281**, 22545 (2006) • [10] Nod1 detects a unique murine peptide from gram-negative bacterial peptidoglycan: S.E. Girardin, et al.; *Science* **300**, 1584 (2003) • [11] Toll-like receptor 2-dependent bacterial sensing does not occur via peptidoglycan recognition: L. H. Travassos, et al.; *EMBO Rep.* **5**, 1000 (2004)

Overview on LPS-forms – “R versus S”

Colony morphology is indicative of the O-glycosylation status. Wild-type bacteria form smooth colonies, synthesize “smooth” or S-form LPS that contain O-polysaccharide chains. S-form LPS also contain R-form LPS molecules in variable proportion depending on culture conditions. So-called “rough-mutants” of gram-negative bacteria synthesize “rough” or R-form LPS. These R-forms lack the O-polysaccharide chains and the core saccharide chain may be present in different stages of completion, giving rise to defined R-classes (e.g. Ra, Rb, Rc, Rd, Re). All R-form LPS are devoid of any S-form LPS [1].

The S-form LPS are commonly the preferred choice for whole animal studies, activating TLR4-positive cells strictly dependent on the presence of membrane-anchored or soluble CD14. They also activate the TLR4/MyD88-independent pathway (TRAM/TRIF: type I IFN) and are therefore selective for classical APC expressing CD14 (e.g. monocytes, macrophages, DC) *in vivo*. The *in vivo* activity of S-form LPS may be modulated by increased levels of LBP (LPS-binding protein) during inflammation [2].

The R-form LPS and Lipid A are primarily used in cellular *in vitro* activation studies. They activate TLR4-positive cells independent of the presence of membrane-anchored or soluble CD14. They do not activate the TLR4/MyD88-independent pathway (TRAM/TRIF: type I

IFN), but they also activate non-classical APC (PMN/mast cells) *in vivo*. They are very useful for *in vitro* cellular activation assays, where CD14/LBP may be absent or only available in limited amounts. A nontoxic nonpyrogenic derivative of Lipid A is Monophosphoryl Lipid A (MPL-A) exhibiting adjuvant properties that may be used in vaccine development [3].

LIT: [1] R-form LPS, the master key to the activation of TLR4/MD-2-positive cells: M. Huber, et al.; Eur. J. Immunol. 36, 701 (2006) ▪ [2] CD14 is required for MyD88-independent LPS signaling: Z. Jiang, et al.; Nat. Immunol. 6, 565 (2005) ▪ [3] Role of innate immune factors in the adjuvant activity of monophosphoryl lipid A: M. Martin, et al.; Infect. Immun. 71, 2498 (2003)

Technical Note

LPS are amphipatic molecules whose hydrophobicity increases with decreasing length of the polysaccharide chain. Therefore Re class LPS and Lipid A are more hydrophobic than Ra class or S-form LPS. The use of Ca²⁺ and Mg²⁺-free buffers is recommended for the preparation of diluted solutions in order to avoid precipitation and loss of activity for Re-LPS, Lipid A and MLP-A. For *in vivo* applications prepare solutions in glucose instead of PBS.




FIGURE 2: Schematic representation of the different LPS chemotypes:
 GlcNac = N-Acetylglucosamine; Glc = Glucose; Gal = Galactose; Hep = Heptose; KDO = 2-Keto-3-desoxyoctonate.
 Adapted from M. Huber, et al.; Eur. J. Immunol. 36, 701 (2006) [1].

In vivo & in vitro Biological Properties of Salmonella minnesota & E. coli LPS

Salmonella minnesota

Product No.	LPS-form	Constituents	Mac	Dend	B Cells	PMN	Mast
ALX-581-009	S-form	O-polysaccharides-Core-Lipid A	++	++	++	+/-	+/-
ALX-581-011	S-form	O-polysaccharides-Core-Lipid A	++	++	++	+/-	+/-
ALX-581-016	R-form (Ra)	GlcNac-Glc-Gal-(Glc-Gal)-(Hep) _{1,2} -Hep-(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-015	R-form (Rb)	Glc-Gal-(Glc-Gal)-(Hep) _{1,2} -Hep-(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-017	R-form (Rc)	Glc-(Hep) _{1,2} -Hep-(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-018	R-form (Rd)	Hep-(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-008	R-form (Re)	(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-201	-	Lipid A	++	++	+++	+++	+++
ALX-581-202	-	Monophosphoryl Lipid A (MPL-A)	++	++	+++	+++	+++

E. coli

Product No.	LPS-form	Constituents	Mac	Dend	B Cells	PMN	Mast
ALX-581-012	S-form	O-polysaccharides-Core-Lipid A	++	++	++	+/-	+/-
ALX-581-013	S-form	O-polysaccharides-Core-Lipid A	++	++	++	+/-	+/-
ALX-581-010	R-form (Ra)	Complete E.coli-core (type II)-Lipid A	++	++	+++	+++	+++
ALX-581-014	R-form (Rc)	Glc-(Hep) _{1,2} -Hep-(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-007	R-form (Re)	(KDO) _{2,3} -Lipid A	++	++	+++	+++	+++
ALX-581-200	-	Lipid A	++	++	+++	+++	+++

CHART LEGEND: Mac = Macrophages, Monocytes ▪ Dend = Dendritic Cells (monocyte-derived) ▪ PMN = Polymorphonuclear Leukocytes (PMN) ▪ B Cells (mouse only) ▪ Mast = Mast Cells
 ▪ +/- = Weak/Absent activation ▪ ++ = Strong activation ▪ +++ = Very strong activation

CpG ODNs & iODNs – Ligands of Toll-like Receptor

- **Activity tested**
- **Potent (0.5-5 µg/ml), endotoxin-free, and selective activators for TLR9 as confirmed with TLR9^{-/-} macrophages and splenocytes**
- **Ultrapure**
- **Easy to handle**

Toll-like receptors (TLRs) are widely expressed recognition receptors of the innate immune system. Four out of ten TLRs identified in humans today recognize nucleic acids, which demonstrates the fundamental importance of microbial DNA and RNA in response to pathogenic microorganisms.

Toll-like receptor 9 (TLR9) recognizes unmethylated CpG motifs in viral and bacterial DNA. Signalling through TLR9 can be accomplished

iODNs – Potent Inhibitors of TLR9 Signalling




FIGURE 3: TLRgrade™ CpG ODN 1668 (Prod. No. ALX-746-001) was added at 50pM/ml to macrophages in a 96-well plate simultaneously together with TLRgrade™ iODN 2088 (Prod. No. ALX-746-250) at the indicated molar excess, cell supernatants harvested after 24 hours and IL-12-p40 analyzed by cytokine ELISA. Similar inhibition rates of >90 (3:1 molar excess) or >95% (10:1) were observed with TLRgrade™ iODN 2088 and CpG ODN 1668, when analyzed for IL-6.

with high efficiency with small, synthetic oligodeoxynucleotides (ODNs) containing optimized CpG motifs, which show great promise as potent vaccine adjuvants and inhibitors of Th2-mediated allergic responses.

Yet uncontrolled TLR9 activation can have deleterious consequences, exacerbating inflammatory tissue damage and increasing sensitivity to toxic shock. Inappropriate TLR9-signalling may also be associated with the promotion of autoimmune diseases like B cell hyperreactivity and anti-DNA antibody production in lupus.

Natural and synthetic DNA sequences have been identified that are able to inhibit immune activation through TLR9. These sequences are derived from diverse sources, including viral sequences, mutated CpG sequences, and repeats of the TTAGGG motif present in mammalian telomeres.

Inhibitory oligodeoxynucleotides (iODNs) composed of TTAGGG multimers reproduce this suppressive activity and block the colocalization of CpG DNA with TLR9 within endosomal vesicles. The mechanism of action is not known, but evidence suggests that inhibitory ODNs do not interfere with cellular uptake of CpG ODNs, nor do they simply compete for binding of TLR9. It was shown that iODNs can affect Th1 priming through the inhibition of STAT1, -3, and -4.

Specific inhibitors of TLR9 would represent a valuable tool for understanding TLR9-mediated responses and might serve as potential therapeutics for some autoimmune diseases.

- **Sterile and pyrogen-free water (1 vial) is now included with every order of CpG ODNs, ODN Controls, and iODNs!**
- **Sterile and pyrogen-free PBS is also available (to be ordered separately).**

Wide Panel of ODNs

Stimulatory ODNs (CpG ODNs)

Three types of stimulatory ODNs have been described, they differ in their sequences and in the type and magnitude of immune responses induced:

- **Type A** are characterized by poly-G motifs with phosphorothioate (PS) linkages at the 5' and 3' ends and a phosphodiester (PO) palindromic CpG-containing sequence in the ODN center. They are very strong inducers of interferon- α (IFN- α) by plasmacytoid dendritic cells (pDC) and especially potent NK cell activators.
- **Type B** are characterized by a full phosphorothioate (PS) backbone with one or more CpG motifs without poly-G motifs. They are weaker inducers of IFN- α , but are potent activators of B cells.
- **Type C** are characterized by a complete phosphorothioate (PS) backbone without poly-G motifs, but also contain palindromic sequences combined with stimulatory CpG motifs.

They offer combined features, i.e. induction of IFN- α by pDC and activation of B cells.

Optimal sequences in ODNs responsible for activating TLR9 vary among species. Recent studies suggest that the species-selectivity attributed to some CpG motifs compared to GpC motifs in control ODNs may only be observed with phosphorothioate bond linkage (Type B and Type C) rather than ODNs containing „natural“ phosphodiester linkages (Type A). The use of TLRgrade™ reagents shows that relative species selectivity for active CpG ODNs often depends on the concentration used.

Control ODNs (GpC ODNs)

Inactive control compounds for CpG ODNs do not stimulate TLR9. They are composed of same sequence as their stimulatory counterparts, but instead of CpG they contain GpC dinucleotides.

Technical Note

Ultrapure TLRgrade™ ODN preparations are selective activators of TLR9 as confirmed with mouse TLR9^{-/-} macrophages and splenocytes. To guarantee consistent and authentic results, the use of endotoxin-free ddWater or PBS is recommended for solubilization. These products have been subjected to multiple rounds of LPS removal by adsorption with activated charcoal.

LT: The removal of ¹⁴C labeled endotoxin by activated charcoal: A.S. Pegues, et al.; Int. J. Artif. Organs 2, 153 (1979)

PBS (endotoxin-free)

ALX-505-007-LD15 1.5 ml
Sterile and endotoxin-free PBS for use with TLRgrade™ or endotoxin-free grade reagents.

ddWater (endotoxin-free)

ALX-505-008-LD15 1.5 ml
Sterile, double distilled and endotoxin-free water for use with TLRgrade™ or endotoxin-free grade reagents.

Selected Literature References

Inhibitory oligonucleotides specifically block effects of stimulatory CpG oligonucleotides in B cells: L.L. Stunz, et al.; Eur. J. Immunol. **32**, 1212 (2002) • Inhibitory oligonucleotides block the induction of AP-1 transcription factor by stimulatory CpG oligonucleotides in B cells: P. Lenert, et al.; Antisense Nucleic Acid Drug Dev. **13**, 143 (2003) • Suppressive oligodeoxynucleotides inhibit Th1 differentiation by blocking IFN- γ and IL-12-mediated signaling: H. Shirota, et al.; J. Immunol. **173**, 5002 (2004) • Inhibitors of TLR-9 act on multiple cell subsets in mouse and man in vitro and prevent death in vivo from systemic inflammation: O. Duramad, et al.; J. Immunol. **174**, 5193 (2005) • Suppressive oligodeoxynucleotides protect mice from lethal endotoxic shock: H. Shirota, et al.; J. Immunol. **174**, 4579 (2005) • Therapeutic potential of oligonucleotides expressing immunosuppressive TTAGGG motifs: D.M. Klinman, et al.; Ann. NY Acad. Sci. **1058**, 87 (2005)

Selected Literature References

Characterization of three CpG oligodeoxynucleotide classes with distinct immunostimulatory activities: J. Vollmer, et al.; Eur. J. Immunol. **34**, 251 (2004) • Cutting edge: species-specific TLR9-mediated recognition of CpG and non-CpG phosphorothioate-modified oligonucleotides: T.L. Roberts, et al.; J. Immunol. **174**, 605 (2005)

Inhibitory/Suppressive ODNs (iODNs)

Inhibitory or suppressive ODNs have been described to neutralize the stimulatory effect of CpG ODNs when added in excess. Two classes of inhibitory ODNs (iODNs) have been described:

- **Class One** is characterized by the presence of methylated CG or unmethylated GC sequences that selectively block CpG-induced immune activation.
- **Class Two** of iODNs contains repetitive TTAGGG motifs patterned after those present in mammalian telomeres.

Stimulatory CpG ODNs & Control GpC ODNs

Now available: **BULK** for in vivo studies!

Product/Sequence	Type	Recommended Species	Prod. No.	Size
ODN 1668 (TLRgrade™) (synthetic) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-001-T100	100 tests
ODN 1668 (TLRgrade™) (synthetic) (BULK) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-051-M001	1 mg
ODN 1720 (TLRgrade™) (synthetic) (Control) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-200-T100	100 tests
ODN 1826 (TLRgrade™) (synthetic) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-002-T100	100 tests
ODN 1826 (TLRgrade™) (synthetic) (BULK) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-052-M001	1 mg
ODN 1982 (TLRgrade™) (synthetic) (Control) 5'-tccatgacgttcctgatgct-3'	Type B	Mouse	ALX-746-201-T100	100 tests
ODN 1585 (TLRgrade™) (synthetic) 5'-ggGGTCAACGTTGAgggggg-3'	Type A	Mouse	ALX-746-003-T100	100 tests
ODN 2118 (TLRgrade™) (synthetic) (Control) 5'-ggGGTCAAGCTTGAgggggg-3'	Type A	Mouse	ALX-746-203-T100	100 tests
ODN M362 (TLRgrade™) (synthetic) 5'-tcgtcgttcgaacgacgttgat-3'	Type C	Human/Mouse	ALX-746-004-T100	100 tests
ODN M383 (TLRgrade™) (synthetic) (Control) 5'-tgctgcttcgaacgacgttgat-3'	Type C	Human/Mouse	ALX-746-204-T100	100 tests
ODN 2216 (TLRgrade™) (synthetic) 5'-ggGGGACGATCGTCgggggg-3'	Type A	Human	ALX-746-005-T100	100 tests
ODN 2216 (TLRgrade™) (synthetic) (BULK) 5'-ggGGGACGATCGTCgggggg-3'	Type A	Human	ALX-746-055-M001	1 mg
ODN 2243 (TLRgrade™) (synthetic) (Control) 5'-ggGGGACGATGCTGgggggg-3'	Type A	Human	ALX-746-205-T100	100 tests
ODN 2006 (TLRgrade™) (synthetic) 5'-tcgtcgtttgtcgtttgtcgtt-3'	Type B	Human/Mouse	ALX-746-006-T100	100 tests
ODN 2006 (TLRgrade™) (synthetic) (BULK) 5'-tcgtcgtttgtcgtttgtcgtt-3'	Type B	Human/Mouse	ALX-746-056-M001	1 mg
ODN 2137 (TLRgrade™) (synthetic) (Control) 5'-tgctgctttgtcgtttgtcgtt-3'	Type B	Human/Mouse	ALX-746-206-T100	100 tests
ODN 2395 (TLRgrade™) (synthetic) 5'-tcgtcgtttcggcgcgccg-3'	Type C	Human/Mouse	ALX-746-020-T100	100 tests

NOTE: Lower case letters indicate phosphorothioate linkage.

Inhibitory ODNs (iODNs)

Product/Sequence	Type	Recommended Species	Prod. No.	Size
iODN 2088 (TLRgrade™) (synthetic) 5'-tcctggcggggaagt-3'	Class One	Mouse	ALX-746-250-T050	50 tests
iODN (ttagg)₄ (TLRgrade™) (synthetic) 5'-tttaggtaggtaggtagg-3'	Class Two	Human	ALX-746-251-T050	50 tests

NEW AT-ODNs see Backcover.

Latest Additions

NEW AT-ODNs!

AT-ODN-1 (endotoxin-free) (synthetic)

ALX-746-024-C100 100 µg

Synthetic. IFN- γ inducing non-CpG ODN of the AT-type, found in the Malaria genome. TLR9-dependent immune activation. **SEQUENCE:** TATAATTTTAATTCCAA-GA. Nucleotides depicted in italics show the corresponding AT-ODN sequence. Includes 1 vial of ddWater (endotoxin-free) (Prod. No. ALX-505-008).

AT-ODN-2 (endotoxin-free) (synthetic)

ALX-746-025-C100 100 µg

Synthetic. *Lactobacillus gasserii*-derived non-CpG ODN of the AT-type. TLR9-dependent immune activation. **SEQUENCE:** TATAATTTTACCAACTAGC. Nucleotides depicted in italics show the corresponding AT-ODN sequence. Includes 1 vial of ddWater (endotoxin-free) (Prod. No. ALX-505-008).

LIT: AT oligonucleotides inducing B lymphocyte activation exist in probiotic *Lactobacillus gasserii*: H. Kitazawa, et al.; Int. J. Food Microbiol. **65**, 149 (2001) • Augmentation of T(H)-1 type response by immunoactive AT oligonucleotide from lactic acid bacteria via Toll-like receptor 9 signaling: T. Shimosato, et al.; BBRC **326**, 782 (2005) • Strong immunostimulatory activity of AT-oligo-deoxynucleotide requires a six-base loop with a self-stabilized 5'-C...G-3' stem structure: T. Shimosato, et al.; Cell. Microbiol. **8**, 485 (2006)

AT-ODN-3 (endotoxin-free) (synthetic)

ALX-746-026-C100 100 µg

Synthetic. IFN- γ inducing non-CpG ODN of the AT-type, found in the Malaria genome. TLR9-dependent immune activation. **SEQUENCE:** TTAACAATTTTAC-CCAAGA Nucleotides depicted in italics show the corresponding AT-ODN sequence. Includes 1 vial of ddWater (endotoxin-free) (Prod. No. ALX-505-008).

LIT: see Product No. ALX-746-025.

NEW Ligands for Intracellular DNA Sensors!

Poly(deoxyadenosine:deoxythymidine) (endotoxin-free) (synthetic)

[Poly(dA:dT) (endotoxin-free) (synthetic)]
ALX-746-022-C050 50 µg

Synthetic. Strong IFN inducer, **independent** of IRF3.

Oligo(deoxyadenosine:deoxythymidine) (endotoxin-free) (synthetic)

[Oligo(dA:dT) (endotoxin-free) (synthetic)]
ALX-746-023-C050 50 µg

Synthetic. Strong IFN inducer, **dependent** of IRF3.

Coming Soon

NEW Endogenous Danger Signal Molecules for TLR4!

Hyaluronic Acid Fragments (<1,500 Da) (TLRgrade™) (Ready-to-Use)

[sHA (<1,500 Da) (TLRgrade™) (Ready-to-Use); Hyaluronan Oligosaccharide (<1,500 Da) (TLRgrade™) (Ready-to-Use)]
ALX-580-004-C250 250 µg INQUIRE!
Isolated from bacteria. **PURITY:** \geq 95%; Activity and endotoxin tested: TLRgrade™. The purity and size range of HA oligomer mixture has been confirmed by HPLC analysis and mass spectrometry. **ENDOTOXIN CONTENT:** <0.02EU/µg. **QUANTITY:** Sufficient for ~100 tests in 100µl assays.


LIT: Oligosaccharides of Hyaluronan activate dendritic cells via toll-like receptor 4: C. Termeer, et al.; J. Exp. Med. **195**, 99 (2002)

Hyaluronic Acid Fragments (1,500-3,000 Da) (TLRgrade™) (Ready-to-Use)

[sHA (1,500-3,000 Da) (TLRgrade™) (Ready-to-Use); Hyaluronan Oligosaccharide (1,500-3,000 Da) (TLRgrade™) (Ready-to-Use)]
ALX-580-005-C250 250 µg INQUIRE!
Isolated from bacteria. **PURITY:** \geq 95%; Activity and endotoxin tested: TLRgrade™. The purity and size range of HA oligomer mixture has been confirmed by HPLC analysis and mass spectrometry. **ENDOTOXIN CONTENT:** <0.02EU/µg. **QUANTITY:** Sufficient for ~100 tests in 100µl assays.

LIT: Oligosaccharides of Hyaluronan activate dendritic cells via toll-like receptor 4: C. Termeer, et al.; J. Exp. Med. **195**, 99 (2002)

Order Your FREE Copy of the
NEW NLR Flyer today at
www.alexis-biochemicals.com!



SWITZERLAND/REST OF EUROPE

ALEXIS CORPORATION
T +41 61 926 89 89
F +41 61 926 89 79
E alexis-ch@alexis-corp.com

NORTH AMERICA

AXXORA, LLC
T (858) 658-0065/1-800-900-0065
F (858) 550-8825/1-800-550-8825
E axxora-usa@axxora.com

BENELUX

10P's AXXORA BVBA
T +32 (0) 3 466 04 20
F +32 (0) 3 466 04 29
E info@10Ps.com

GERMANY

AXXORA DEUTSCHLAND GmbH
T 07621 5500 522
Toll Free 0800 253 94 72
F 07621 5500 523
E axxora-de@axxora.com

UK

AXXORA (UK) Ltd.
T 01949 836111
F 01949 836222
E axxora-uk@axxora.com

International Distributors: Australia Sapphire Bioscience (02) 9698 2022 Austria Eubio (01) 8950145 Bangladesh Future Business Vision (02) 863 1173 Bosnia & Herzegovina A-Z Consulting +386 1 433 63 22 / +386 1 230 18 84 Brazil Biogenity (011) 3666 3565 / Sellex (011) 5506 4646 Canada Cedarlane Laboratories (289) 288-0001/1-800-268-5058 Chile Biocant (2) 683 2437 China Beijing Bitab Biotech (010) 8201 5225 / Boppard 21 6288 4751 (Shanghai) 20 8732 6381 (Guangzhou) / ITS China (021) 6481 4428/98 Colombia Lab Scientific +1 305 716 9922 Czech Republic Genetica (02) 7270 1055 Denmark Medinova Scientific 3956 2000 Ecuador Lab Scientific +1 305 716 9922 Egypt New Test For Scientific Service 03-358-3543 Estonia In vitro Eesti 630 65 20 Finland Nuppulinan Laboratoriopalvelu (02) 07920350 France CovalAb 0437 654 236 / Coger (01) 45 33 67 17 Greece SB Bio-technology Suppliers (210) 823 3373 Hong Kong Boppard 2799 9019 Hungary Biomarker 28 419 986 India Hysel India 011-2622 78 01/02/03/04 / Imgenex India (01674 274 3265/(01674 329 6544 Indonesia ITS Indonesia (021) 451 6222 Iran Hormoz Pajohan Lab. Equipment (021) 888 3444 Iraq Alwojoh (01) 515 8350 Ireland Alpha Technologies 045 865 440 / 0149 62 422 Israel Almog Diagnostic (03) 977 3390 Italy Vinci-Biochem 0571 568147 Japan BioLinks 03 5443 6891 Korea Chun Yang Tech (02) 929 8071 Lithuania & Latvia In vitro Eesti +372 630 65 20 Malaysia Interscience (03) 57 40 9888 Mexico Consultoria de Laboratorios (0155 1163 8840 The Netherlands 10P's Axxora 076 5425 184 New Zealand Sapphire Bioscience +61 2 9698 2022 Northern Ireland Alpha Technologies 028 28260558 Norway AH Diagnostics (23) 23 32 60 Pakistan The Worldwide Scientific (0142-755-2355 Poland Biomibo (022) 872 0797 Portugal Baptista Marques (21) 722 06 60 Romania Medist (21) 411 5003 Russia Chimmed 095 728 4192 Singapore ITS Science & Medical (06) 273 0898 Slovenia A-Z Consulting (01) 433 63 22 / (01) 230 18 84 South Africa Southern Cross Biotechnology (021) 671 51 66 Spain Grupo Taper 916 596 520 Sweden Immunkemi 08-58361500 / In vitro 08-306010 Syria New-Med Technology (11) 8827 1717 Taiwan Cashmere Scientific Company 0800 222 095/02-2567 5682 Thailand ITS Thailand (02) 308 0611 / Theera Trading (02) 412 5672 / (02) 418 1068 / S.M. Chemical Supplies (02) 542 1037 Turkey Tokra (312) 395 6009 Uruguay Lab Scientific +1 305 716 9922 Venezuela Lab Scientific +1 305 716 9922 Vietnam ITS Vietnam (08) 9255 232

ALEXIS[®]
BIOCHEMICALS

Available on

www.axxora.com